




SCALE 1:1

SECTION A-A



ALL PRINTED COPIES ARE FOR REFERENCE ONLY

00	20/01/2023	OUB	Creation			
Rev.	Date	Designer	Modification			
R&D Approv.: JJA Date: 20/01/2023			QA Approv.: Date:		Cust. Approv.: Date:	Norms: ISO 594-1/2
Material: ABS					Color: Clear	
Material code: M0071					Revision: 00	
Reference: AAL01050 Designation: MALE LUER LOCK / FEMALE LUER LOCK – CLEAR ABS					General tolerances: .X ±0.3 .XX ±0.13 Units: mm ° grs	
					Scale: 4:1	
					 PROMEPLA	
					FOR035-4 P. 1 / 1 A4	