




SCALE 1:1

COUPE A-A



00	21/02/2022	OUB	Creation			14046
Rev.	Date	Designer	Modification			N° FOR073 ou FOR021
Verif R&D : PPE Date : 22/02/2022		Verif QA : CGI Date : 25/02/2022		Verif client : N/A Date :	Norms : N/A	
Material : PP			Color : Natural			Revision : 00
Material code : M0598						State : VALID
Reference : ABF16050 Designation : FEMALE LUER LOCK FOR TUBE ID 1,6MM				General tolerancies : N/A		Scale : 5:1
				Units : mm, °, grs		 FOR035-3 P. 1 / 1 A4