



ORIGINAL

**NOTA :**

Rep	Matière	Code matière
1	PMMA	M0115
2	Gortex 0.03µPTFE+PES 0.2µ	M0215 + M0249

Product specifications	
Priming volume	2 ml (Max)
Minimum water flow rate	20 ml/min @ 36 inch
Minimum water bubble point	40 psi mini for 10 sec
Effective filtration area	11 cm <sup>2</sup>
Maximum operating pressure	45 psi (Max)

Rev.	Date	Designer	Modifications	Internal Validation	Quality Validation
2	30/07/2012	BOUJMADA	Mise à jour matière ; M0107 devient M0115, M0332 devient M0215, cotation, design et mise au format SW		
1	24/03/2010	RANDAZZO	Remplacement du certificat matière M0330 par M0332		
0	11/04/2006	BOUCHARD	Creation		

Dessinateur : BOUCHARD Date : 11/04/2006	Visa : <i>Fred Dauchin</i> Date : <i>16 08 12</i>	PROMEPLA	Echelle : 1:1 Unités : mm, °, grs	Tolérances générales :
Matière : VOIR NOTA			Normes :	
Code matière : VOIR NOTA			P. 1 / 1	
Référence : EDB04990 Désignation : FILTRE ADULT I.V. HYDROPHILE 0.2µ POUR TUBE 2xØ1 3mm			Révision : 2	 FOR035-0