

AM 22 W17
EVICOM AM/22-W17

Certificate Number : M0240 dated:12/12/2000
printed:27/03/2015

Technical data

Specific weight	ISO 1183	1.31 g/cm ³
Melt flow rate (MFR 190°C/ 20.325kg)	ASTMD 1238	28 g/ 10min
Hardness (Shore D)	ISO 868	80 SHD
Tensile strength at yield	ISO 527	34 N/mm ²
Elongation at break	ISO 527	230 %
Vicat softening point (49N)	ISO 306	54 °C

Main Medical Compliance:

ISO 10993-3/4/-5/-10/-11
USP XXVIII/ Class VI

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: PVC Rigid- LATEX Free- NO BPA- NO ESB RISK- No Phthalates- BPB and PCB free- REACH compliance (Directive 1907/2006/EC) - SVHC list free (03/2009)- RoHS Compliance (Directive 2002/95/EC)-MSDS- other ref: 9700200000