



PVC VJ 375/1 CLEAR42

Certificate Number : M0258 dated:05/02/2001
printed:27/06/2017

Technical data

EXTRACT CAD'M LEVEL	CTM 59	0.000 PPM
HEAT STABILITY 180°C	CTM 4	195 MINS
MACKLOW FINISH	CTM 6	4.5
MACKLOW PRESSURE	CTM 6	90 KG/CM2
NMR	CTM 14	4.89
EXTRACT LEAD LEVEL	CTM 59	0.000 PPM
SPECIFIC GRAVITY	CTM 9	1.210
EXTRACT TIN LEVEL	CTM 59	0.000 PPM
UV ABSORPTION	CTM 67	0.000
EXTRACT ZINC LEVEL	CTM 59	1.103 PPM

Main Medical Compliance:

E.P. VERSION VI 1.2.1.1
USP CLASS VI

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: PVC SOUPLE-presence of DEHP (31%)-No BPA