



M.G.S.

## HS91 CRISTALLO

Certificate Number : M0279 dated:02/02/1999  
printed:27/06/2017

### Technical data

Specific weight		1.26
Load at break		200 MPa
Elongation at break		270 %
Hardening temperature		-5 °C
Acidity(NaOH N/100)/100 ml		0.30 ml

Main Medical Compliance:

European Ph. III Ed.

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization \*: *\*must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: **PVC Souple Shore A 91**