

NAKAN FEM671N T08

Certificate Number : M0685 dated:31/10/2012
printed:25/07/2017

Technical data

Density	ISO 1183	1.21 kg/dm ³ ±0.02
Hardness (Shore A)	ISO 868	67 SHA ±2
Tensile strength at break	ISO 527	12- 16 MPa
Elongation at break	ISO 527	370- 410%

Main Medical Compliance:

USP Class VI
ISO10993-4/ -5
European Pharmacopoeia 3.1.1.1

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: Soft PVC medical grade- Plastizicer DOP (2-ETHYLHEXYL)
(concentration: 5-50%)- SVHC list (19/11/2012)- MSDS