

NAKAN ATA 43

Certificate Number : M0713 dated:18/05/2010
printed:26/07/2017

Technical data

Density	ISO 1183	1,34 kg/dm ³ ±0.02
Hardness (Shore D)	ISO 868	81 ±2 SHD
Tensile strength	ISO 527	30 - 34 MPa
Elongation at break	ISO 527	35 - 75 %
Vicat point	ISO 306	65°C ±2

Main Medical Compliance:

ISO 10993-4/ -5/ -10
USP Class VI
USP XXVI/ XXVII/ XXVIII
ASTM F750-82

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: PVC for medical application- LATEX Free- No BPA- No Phthalates-
No ESB RISK- RoHS compliance- REACH compliance (Directive
1907/2006/EC)- MSDS- Other ref.: NAKAN ATA 43 GMS C.065