

TERLUX 2822HDCertificate Number : M0847 dated:17/10/2011
printed:14/11/2014**Technical data**

Density	ISO 1183	1080 kg/m ³
Melt volume rate (MVR)(220°C/ 10kg)	ISO 1133	2 cm ³ / 10min
Tensile stress at yield	ISO 527	48 MPa
Tensile strain at yield	ISO 527	4 %

Main Medical Compliance:USP Class VI
ISO 10993-5
European Pharmacopeia 3.2.2
DMF N°25140 (FDA)**Recommended Uses:**

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: MABS- LATEX Free- No BPA- NO ESB RISK- No Phthalates- SVHC list free (16/06/2014)- REACH compliance (Directive 1907/2006/EC)- RoHS Compliance (Directive 2011/65/EC)- MSDS