

Certificate Number : M1195 dated:29/01/2015
printed:15/08/2017

Technical data

Density	ISO 1183	1.25 kg/ dm3 ±0.02
Hardness (Shore A)	ISO 868	93 SHA ±2
Tensile strength at break	ISO 527	18-22 MPa
Elongation at break	ISO 527	280-320%

Main Medical Compliance:

USP Class VI/ USP XXVIII
ISO10993-4/ -5/ -10

Recommended Uses:

- Extrusion
- Moulding
- Blowing

Granule color:

- Natural
- Red
- Blue
- White

Other:

Recommended Sterilization *: **must be validated by the customer with his sterilization cycle*

- E.O.
- Gamma radiation
- Steam

Note: Soft PVC DEHP Free- Plastizier TEHTM (24%)- LATEX Free- No BPA- No ESB RISK- No Phthalates- RoHS compliance (Directive 2011/65/EU)- SVHC list free (05/04/2013)- REACH compliance (Directive 1907/2006/EC)- MSDS