

NOTA :

Rep	Matière	Code matière
1	PC	M0260
2	PC	M0260
3	Silicone	M1001
	Lubrifiant	M1000

PRODUCT SPECIFICATION	
Flow rate	6 l/h at 1 meter head height
Priming volume	0.04 ml
Back pressure	60 PSI (tested with water)
Negative pressure	-2.9 PSI
Displacement	-0.008 < > 0
Sterilization	ETO and Gamma radiation



APPROBATION / APPROVAL		
Client / Customer :	Promepla Resp. BE :	Promepla Qualité :
	<i>[Signature]</i>	
11 08 15		
Nom / Name, date & signature		


1	09/03/2015	BOUJMADA	Mise à jour désignation	FOR021 : 6686
0	11/12/2014	RANDAZZO	Creation (Proto 111214C)	FOR073 : 8359
Rev.	Date	Designer	Modifications	N° Fiche FOR073 ou FOR021

Normes : ISO 594-2

Matière : VOIR NOTA

P. 1 / 1

Echelle : 2:1



Code matière : VOIR NOTA

Revision : 1

Format : A4

Reference : ICB11990

Designation : SITE CLOS A ACTIVATION LUER - LLF/ LLM

Tolérances générales :

Unités : mm, °, grs

PROMEPLA

FOR035-1